

EMERGENCY

Emergency bag



We have concentrated on the functionality of the product: a lightweight bag, easy to handle and easy to clean. Eliminating the superfluous, all you have left is the handiness of the product and the convenience of its price. Made of nylon, it has a central compartment where you will find on the inside a pre-modelled and semi-rigid compartment with three dividers to store a small oxygen tank and other resuscitation devices. The closure of the top flap is fast and safe thanks to the Velcro® strap all around the front border. There are also two side pockets with zipper.

The Emergency stands out for the global character of its range: it is available in fact in different languages.

TECHNICAL DATA

Material: Ristop 600 D
 Dimensions: 525 x 215 x h230 mm
 Weight empty: 700 g
 Weight complete: 6 kg
 Capacity: 26 L
 Color: Orange / Black

Content

COMPLETE EMERGENCY BAG "UE"
 ADULT PVC RESUSCITATOR W/POP OFF
 SPENCER GUARANTEE 1 COLOR
 ORANGE TUFF CUT UNIVERSAL SCISSORS
 DGX 410 ANEROID SPHYGMOMANOMETER W/ADULT BLUE CUFF
 DGX 110 BLACK BASIC PHONENDOSCOPE - NO LATEX
 PLASTIC BLUE PEN LIGHT W/2 BATT. LOGO AND BLISTER
 WHITE ADHESIVE LABEL 60x60mm (2000 pcs/roll)
 MAGILL ADULT FORCEPS 268 mm.
 TUBULAR TOURNIQUET IN RUBBER LATEX "CE" 350mm
 BOX E AVANA cm.43x34x28h 3 WINGS
 O2 BOTTLE 0,5lt UNI W/VALVE
 ANTI-CRUSH TUBE FOR OXYGEN
 PISTON PRESSURE REDUCER B3,5-6l/m-SHORT TANG-W/
 COVER GAUGE UNI CONNECTION
 DISPOSABLE NON-STERILE GUEDEL mm.60
 DISPOSABLE NON-STERILE GUEDEL mm.70
 DISPOSABLE NON-STERILE GUEDEL mm.80
 DISPOSABLE NON-STERILE GUEDEL mm.90
 CUFFED MASK SIZE 3 PVC COLOUR CODED WITHOUT VALVE
 CUFFED MASK SIZE 4 PVC COLOUR CODED WITHOUT VALVE
 STAINLESS STEEL YOUNG TONGUE FORCEPS
 JAW SPREADER MOUNDER
 RES Q VAC MANUAL ORANGE SUCTION PUMP
 EMERGENCY BAG W/"EMERGENCY" LOGO
 DEPLIANT BAGS CROSS SELLING
 ISOTHERMAL SHEET ADULT cm.160x220 W/LOGO

Emergency CB08000