

# ACTAR D-fib®



- From the company who pioneered the One Manikin Per Student Concept in CPR Training.
- Compatible with all AED electrodes for use with AED Trainer.
- Jaw thrust opens the airway without head tilt to simulate spinal injury patient.
- Realistic resistance for proper chest compressions and chest rise.
- Incorporates all anatomical landmarks needed to teach CPR.
- Patent pending design permits stacking of manikins for easy transportation and storage.
- Lung/mouth shield eliminates the need to soak manikin parts.

Product Description		Actar D-fib®	
Ordering Information	Product Code	Description	Qty Per Case
	AC01	ACTAR D-fib® Compact 1™ Includes: 1 ACTAR D-Fib manikin, 5 disposable lungs, back pack	1
	AC05	ACTAR D-fib® Compact 5™ Includes: 5 ACTAR D-fib manikins, 50 disposable airway lungs, nylon back pack	5
	AC10	ACTAR D-fib® Compact 10™ Includes: 10 ACTAR D-fib manikins, 50 disposable lungs, carry bag	10
	AC25	ACTAR D-fib® Compact 25™ Includes: 25 ACTAR D-fib manikins, 100 disposable lungs, rolling cart with nylon cover	25

www.vital-signs.com  
 20 Campus Road, Totowa, NJ 07512, USA • TEL: 800-932-0760 • FAX: +1-973-790-3307  
 Authorized European Representative: Vital Signs, Ltd., 13-14 Eldon Way, Lineside Industrial Estate, Littlehampton, West Sussex, UK



vital signs