

Robust band for the stabilization of head


Rugged handles covered in Santoprene® to improve the grip and guarantee comfort for the operator

Adjustable belt with traditional buckle to give assurance to the patient no matter what their build

Double stitching to ensure durability

Compact woven cotton canvas for heavy loads

Seamless supporting cord to guarantee secure raising and lowering



Available in the following colours



NRT Neil Robertson type

A rescue stretcher specifically designed for the recovery of victims in difficult and cramped situations, especially where the patient has to be raised and carried to safety. Made with semi-rigid splints and high resistant cotton. It guarantees utmost safety for the patient also during lifting.

Thanks to a robust construction and a practical storage system it becomes indispensable for rescue in:

- mines
- industries
- building sites
- ships

Dimensions: 1150 x 20 x h1900 mm

Weight: 8 kg

Load capacity: 172 kg

NRT white

ST30150A

NRT yellow

ST30152A

STX 51 transport bag

ST30151B

S	CE			LXF	
E	HP	B	150 x 100 mm	col 01	min 01
I	interest	FR	MC	MI	MR

