

SPENCER 122

Four fold stretcher



Thanks to reduced clutter, the stretcher is used as an emergency stretcher and also in conditions where storage space lacks.
Complete with two orange quick-release restraints
Reflective logo to increase nocturnal visibility

TECHNICAL DATA

Length open	2010 mm
Length closed	1000 mm
Width open	550 mm
Width closed	110 mm
Height closed	140 mm
Height open	170 mm
Weight	6 kg
Load capacity	170 kg
Handles	4
Support feet	4

ST00122A

SPENCER 122
FOUR FOLD STRETCHER
ORANGE



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>